

## MITSUBISHI GAS CHEMICAL

# MULTILAYER PLASTIC VIALS & SYRINGES FOR BIOLOGICS

Here, Takuya Minezaki, MD, Research Manager, and Tomohiro Suzuki, Associate General Manager, both of Mitsubishi Gas Chemical, provide an overview of the company's OXYCAPT™ Vials and Syringes, explaining how their multi-layered polymer material delivers the advantages of both COP and glass, and is especially suitable for biotherapeutics.

Based on its technologies and experiences, Mitsubishi Gas Chemical (MGC) has successfully developed multilayer plastic vials and syringes called OXYCAPT™ (see Figure 1).

OXYCAPT™ consists of three layers (Figure 2): the drug contact layer and outer layer, both made from cyclo-olefin polymer (COP); and the oxygen barrier layer, which is made from a proprietary, novel polyester. This enables OXYCAPT™ to offer:

- excellent oxygen barrier properties
- high water vapour barrier properties
- excellent UV barrier properties
- very low extractables
- high pH stability
- low protein adsorption and aggregation
- a silicone-oil free barrel
- high transparency

- high break resistance
- easier disposability
- lighter weight.

As regulatory authorities have reported, there are problems with some existing glass and plastic vials and syringes. For example, glass suffers from breakage, delamination, etc, and traditional plastic does not represent a sufficient oxygen and UV barrier. Especially in glass, the US FDA has pointed out such problems and reported more than 50 recalls. To address these problems from glass, many suppliers have launched plastic vials and syringes as alternatives, but their oxygen barrier properties have not always met customer demands.

OXYCAPT™ represents a plastic material with oxygen barrier properties almost the same as glass, and more than 100



**Takuya Minezaki**  
Research Manager  
T: +81 463 21 8627  
E: takuya-minezaki@mgc.co.jp



**Tomohiro Suzuki**  
Associate General Manager  
T: +81 3 3283 4913  
E: tomohiro-suzuki@mgc.co.jp

Figure 1: OXYCAPT™  
Vials and Syringes.



**Mitsubishi Gas Chemical Company, Inc**  
Mitsubishi Building  
5-2 Marunouchi 2  
Chiyoda-ku  
Tokyo 100-8324  
Japan

[www.mgc.co.jp/eng](http://www.mgc.co.jp/eng)




Figure 2: The multilayer structure of OXYCAPT™ Vials and Syringes.




Figure 3: Graph comparing oxygen permeation properties of OXYCAPT™ with COP and glass.




Figure 4: Graph comparing UV transmission properties of OXYCAPT™ with COP and glass.

“OXYCAPT™ vials and syringes are produced by co-injection moulding technology. This technology has been applied to beverage bottles for many years, but we are the first company that has succeeded in developing multilayer plastic syringes.”

times better than COP (Figure 3). According to our internal studies using antibodies, OXYCAPT™ also outperformed both glass and COP in terms of preventing oxidation. Biologics are often vulnerable to oxidation, and OXYCAPT™ can contribute to the stability of such oxygen sensitive drugs.

OXYCAPT™ UV barrier properties also compare very favourably indeed to other materials. For example, about 70% of UV (300 nm) light transmits through glass and COP, whereas only 1.7% transmits through OXYCAPT™ (Figure 4). We have confirmed that this feature of OXYCAPT™ also contributes to stability of biologics.

As a barrier to water vapour, OXYCAPT™ cannot outperform glass. Nonetheless, it is comparable with COP and, like COP, OXYCAPT™ comfortably meets with the ICH guideline for water vapour barrier properties.

Studies have shown extremely low extractables with OXYCAPT™. One study was conducted to confirm volatile, semi-volatile and non-volatile impurities from OXYCAPT™. Water and four solutions (50% ethanol, NaCl, NaOH and H<sub>3</sub>PO<sub>4</sub>) were used and impurities were measured by gas chromatography-mass spectrometry (GC-MS) and liquid chromatography-UV spectroscopy-mass spectrometry (LC-UV-MS) after 70 days at 40°C. Compared with the control, no impurities were detected in any of the OXYCAPT™ containers. A second study was conducted to measure inorganic extractables from OXYCAPT™. The level of extractables was similar to those from COP, which is well-known as an extremely pure polymer, and less than that of Type I glass (Figure 5).

The OXYCAPT™ Syringe consists of tip-cap, barrel, polytetrafluoroethylene (PTFE)-laminated stopper, and plunger rod. Although a very small amount of silicone-oil




Figure 5: Graph comparing inorganic extractables levels from OXYCAPT™ with those from COP and glass.

“Vials were dropped from 150 cm height to confirm break resistance. All of the glass vials were broken, but zero breakage was observed with OXYCAPT™ vials.”

is sprayed on stoppers of the OXYCAPT™ Syringe, no silicone oil is baked-on to the barrel. According to MGC internal studies using antibodies, this feature leads to much less protein aggregation compared with Type 1 glass syringes.

OXYCAPT™ vials and syringes are produced by co-injection moulding technology. This technology has been applied to beverage bottles for many years, but we are the first company that has succeeded in developing multilayer plastic syringes. We

have also successfully developed inspection methods for the oxygen barrier layer such that all containers are 100% inspected by state-of-the-art inspection machinery.

MGC can offer bulk vials, ready-to-use (RTU) vials and RTU syringes. Regarding the RTU products, vials and syringes are provided in ISO based nest & tub formats. The nests & tubs are mainly sterilised by gamma irradiation. MGC offers 2, 6 and 10 mL vials, and 1 mL “Long” and 2.25 mL syringes, and the company is able to provide samples for initial testing free of charge.

Each polymer meets the requirements of USP661, USP87, USP88, EP and has been filed in the FDA’s drug master file (DMF). The vials and syringes are also compliant with each pharmacopoeia and have been filed in the DMF. The syringes are produced and controlled in accordance with ISO 13485.

In recent studies of lyophilisation and cold storage resistance, after being filled with albumin solution, vials were primarily frozen at -50°C for six hours, secondary dried at 4°C for 48 hours, and finally dried at 25°C for nine hours. We measured oxygen barrier, appearance and dimensions, and found OXYCAPT™ maintained its properties and dimensions before and after lyophilisation (Figure 6).


In addition, vials were dropped from 150 cm height to confirm break resistance. All of the glass vials were broken, but zero breakage was observed with OXYCAPT™ vials (Figure 7). Vials were also stored at -80°C to confirm cold storage resistance. When



UNPARALLELED TIGHT  
TOPIC FOCUS, EVERY ISSUE  
FOR FOURTEEN YEARS

[www.ondrugdelivery.com](http://www.ondrugdelivery.com)

## Oxygen Barrier



## Appearance & Dimensions

| 10 mL vial | After Lyophilisation |
|------------|----------------------|
| Appearance | No Change            |
| Dimensions | No Change            |

Figure 6: Resistance to lyophilisation – oxygen barrier properties (top) and dimensions and appearance (bottom) maintained.

frozen OXYCAPT™ vials were dropped from 150 cm height, no breakage was observed. To understand any influences on long-term cold storage, the same studies will be conducted over six-month and two-year periods.

Our targeted application for OXYCAPT™ is biologic therapeutics. As the ICH Guideline “Stability of Biotechnological/Biological Products Q5C” mentions, oxidation is one of the causes of protein instability. The excellent profile of OXYCAPT™ with regard to oxygen and UV barrier properties promises to contribute to stability of biologic products.

In addition, we believe OXYCAPT™ can be applied to epinephrine, because it is well known to be an oxygen sensitive drug. Glass syringes, having problems with breakage, are not ideal for emergency drugs, so some suppliers have tried to develop new pen injectors made of plastic.

Customisability is another feature of plastic, and MGC can customise OXYCAPT™ containers to meet specific requirements.

## CONCLUSION

OXYCAPT™ has been developed to overcome the current problems the pharmaceutical industry is experiencing with syringes and vials made from traditional materials. OXYCAPT™ combines the benefits and features of COP, such as high water vapour barrier properties, high break resistance, very low extractables and low protein adsorption, with the high oxygen & UV barrier properties associated with glass. In the rapidly growing biologics sector, OXYCAPT™ vials and syringes deliver substantial benefits.




Figure 7: Zero OXYCAPT™ breakages during drop test, compared with 100% breakage with glass.

## ABOUT THE COMPANY

Mitsubishi Gas Chemical (MGC) develops and manufactures chemical products ranging from basic to fine chemicals and performance materials. It comprises the Natural Gas Chemicals Company, the Speciality Chemicals Company, the Aromatic Chemicals Company and, of particular relevance to those interested in drug delivery, the Information & Advanced Materials Company and the Advanced Business Development Division.

MGC established its Advanced Business Development Division in 2012 as a centre for continually creating new businesses across the fields of medical/food, information and communications, mobility, energy and infrastructure. The group has developed OXYCAPT™ Vial & Syringe as an alternative to glass containers for parenteral primary packaging.


## ABOUT THE AUTHORS

**Takuya Minezaki** joined Mitsubishi Gas Chemical in 2005. He specialises in the synthesis and injection moulding of polymers, and has worked in the development of oxygen absorbing polymers. Since 2017, he has worked on the development of new pharmaceutical containers by utilising oxygen absorbing resins. He is also in charge of the development of secondary materials, such as plunger stoppers and tip-caps. Mr Minezaki holds a Master of Engineering in Applied Chemistry from the University of Tokyo, Japan.

**Tomohiro Suzuki** joined Mitsubishi Gas Chemical in 1998. He belonged to the Oxygen Absorbers division until 2011, after which he was transferred to Advanced Business Development division in 2012, to be a member of the OXYCAPT™ development team. Since then, he has been in charge of marketing for OXYCAPT™ Vial & Syringe. His current position is Associate General Manager.

# OXYCAPT™ Plastic Vial & Syringe

## Multilayer Structure



Water Vapor Barrier Layer  
(COP)

Oxygen Barrier Layer  
(New Polymer)

Drug Contact & Water Vapor Barrier Layer  
(COP)



- ✓ Excellent Oxygen Barrier
- ✓ High Water Vapor Barrier
- ✓ Low Extractables & High pH Stability
- ✓ High Break Resistance & Lightweight
- ✓ Excellent UV Barrier
- ✓ High Transparency
- ✓ Silicone Oil Free Barrel
- ✓ Low Protein Adsorption & Aggregation
- ✓ Suitable for Biologics
- ✓ Customizable



mitsubishi gas chemical

Mitsubishi Gas Chemical Company, Inc.  
<https://www.mgc.co.jp/eng/products/abd/oxycapt.html>

Mitsubishi Gas Chemical America, Inc.  
<http://www.mgc-a.com>

E-mail: [nb3.pharmapackage@mgc.co.jp](mailto:nb3.pharmapackage@mgc.co.jp)