



OWEN MUMFORD

Pharmaceutical Services


AIDAPTUS®: OWEN MUMFORD INTRODUCES ITS NEXT-GENERATION PLATFORM AUTOINJECTOR

In this article, Michael Earl, Director, Pharmaceutical Services at Owen Mumford, introduces the next-generation Aidaptus autoinjector, explaining how this latest platform device from Owen Mumford can contribute to a faster speed-to-market, increased flexibility during drug development and improved patient outcomes.

Autoinjectors have been used as a convenient and effective means of delivering subcutaneous injections since the 1980s, and recent healthcare trends have led to an increase in their use. Autoinjectors allow patients to administer their own medication outside of acute care facilities, without the need of a healthcare professional. This helps to address a number of healthcare concerns, including an ageing population, the increase in chronic conditions among patients and pressure on healthcare systems and budgets.

A further driver of autoinjector use is the emergence, and subsequent increase, of biologics and biosimilars over the past years. More low volume formulations are being developed for these drugs to allow administration via the subcutaneous route, facilitating patient self-administration. The ensuing increase in demand for autoinjectors has created greater competition in this area, with simplicity and ease of use for the patient being a major differentiator between devices.

Figure 1: Aidaptus – Owen Mumford's next-generation autoinjector.



THE PLATFORM APPROACH

Speed-to-market is a key consideration in the drug development process, meaning that commissioning custom devices for specific drugs may not always be the most cost-effective or efficient solution. In more recent years, pharmaceutical companies have tended to opt for platform devices, which are a quicker and more flexible solution. The base platform device is designed with a wide envelope of possible delivery capabilities so that it can be used for a variety of formulations, reducing the level of testing required at the device selection stage.

However, the challenge with platform devices is that they are not designed with a specific user group in mind; the devices must therefore accommodate the needs of multiple potential patient groups. This broad accommodation is critical to



Michael Earl
Director, Pharmaceutical Services

**Owen Mumford
Pharmaceutical Services**
Brook Hill
Woodstock
Oxfordshire
OX20 1TU
United Kingdom

www.ompharmaservices.com

reassure pharmaceutical businesses that any use-related risk factors have been identified and addressed regardless of the patient demographics. Device developers must therefore devise inclusive testing strategies that cover patients with different levels of physical and cognitive abilities.¹

BRINGING AIDAPTUS TO MARKET

A pioneer in the development and manufacture of autoinjectors, Owen Mumford Pharmaceutical Services has launched a new platform autoinjector offering next-generation benefits of flexibility and versatility (Figure 1). As with all other platform devices in the Owen Mumford Pharmaceutical Services range, the Aidaptus® disposable autoinjector was designed with a focus on ease of use and patient comfort.

However, the covid-19 pandemic began during the latter stages of the Aidaptus' development, meaning that final human factors testing could not proceed as usual. It is only through first-hand user feedback on a range of device prototypes that usability and human factors issues can be properly understood, so, to be able to move forward, the human factors team decided to prepare a new study design that would comply with the restrictions in place. For example, the prototypes were all packed at least three days in advance of the study by engineers wearing PPE, and were kept sealed until needed. Rather than handing out each prototype as usual, the moderator instructed participants on which prototypes were needed at each stage and which actions they needed to take.² This carefully adapted study provided the key data to allow the team to meet human factors requirements and to launch the product this autumn.

DESIGN BENEFITS FOR PATIENTS

Understanding the stages of needle insertion and medication delivery is key to helping patients to manage their injection

“Needle insertion and dose delivery take place as two separate phases, with each action controlled independently. This means that the delivery of the medication does not take place until the needle is fully inserted into the patient’s skin, helping to prevent drug spillage prior to injection.”

process successfully and to be confident that they have administered their drug correctly. For this reason, the Aidaptus autoinjector provides patients with audible and visual notifications during the injection procedure (Figure 2). There are audible clicks at the start and end of dose delivery, as well as a bright yellow plunger rod clearly visible through the viewing window to confirm end of dose. The window is large enough for patients to easily see the drug before administration, and to check the drug’s clarity and colour. There are no other internal mechanisms visible through the window, thereby providing an unobstructed view, as well as confidence, to the user. The needle itself is hidden before, during and after use to prevent needle exposure – a feature that may be particularly beneficial for patients with needle-phobia. On completion of the injection, the safety shroud locks in place to ensure that the needle cannot be exposed and present a risk to the user before disposal.

HELPING TO PREVENT WET INJECTIONS

An issue occasionally experienced while using autoinjectors is the occurrence of “wet injections” – patients may remove the autoinjector from the injection site too early in the drug delivery process, resulting in drug spillage and wastage. Aidaptus’s



Figure 2: Aidaptus provides both audible and visual feedback to give patients confidence in their self-injection.

design mitigates this through automatic pressure-activated needle insertion that creates a consistent user experience for all injections. Needle insertion and dose delivery take place as two separate phases, with each action controlled independently. This means that the delivery of the medication does not take place until the needle is fully inserted into the patient’s skin, helping to prevent drug spillage prior to injection.

This dual-phase action also helps to prevent syringe breakage, which can be caused by the strong springs typically used to deploy the needle for skin insertion and the plunger for dose delivery. Furthermore, before the product even reaches the patient, the self-adjusting plunger helps to prevent breaches of container closure integrity by

“For maximum flexibility, the Aidaptus autoinjector is compatible with either a 1 or 2.25 mL prefilled syringe, with a minimal number of changed parts, whilst maintaining its small, discreet size (162 by 18 mm).”



Aidaptus[®]

2-step single-use auto-injector platform

Versatile design intuitive delivery

Fill Volume Flexibility

Aidaptus[®] readily adapts to a range of different drug fill volumes with no changed parts, using a self-adjusting plunger.



No change parts required for fill volume changes

* In addition to an air bubble & overfill

Want to know more?

Find out more about our new innovative Aidaptus[®] auto-injector by scanning the QR code or visiting www.ompharmaservices.com/odd



 **OWEN MUMFORD**
Pharmaceutical Services

Aidaptus[®] is a registered trademark of Owen Mumford Ltd, ©2021

OMPS/odd/ad/ob/1021/07

Available now