



MICRONISATION SUCCESS STARTS WITH SOLID-STATE UNDERSTANDING

Catalent
Pharma Services™

Dr Derek Beauchamp and John Kusnierz, both of **Catalent**, explain how solid-state characterisation underpins the behaviour of a drug substance under processing conditions, and consider the importance of micronisation in improving bioavailability and drug product performance.

As molecules become more complex and solubility challenges more common, particle engineering, particularly micronisation, has become a core strategy to enhance bioavailability and optimise drug product performance. For BCS Class II compounds, reducing mean particle diameter to the low-micron range increases specific surface area and accelerates dissolution kinetics, consistent with the Noyes–Whitney relationship.^{1,2} In highly potent API programmes, micronisation also enables precise dosing and content uniformity when API loading is minimal.

However, micronisation is not a neutral unit operation. The mechanical energy introduced during size reduction can alter the solid-state form of the API by way of polymorphic transitions, partial amorphisation and surface disorder.^{3,4} These changes may not be immediately visible, but they can directly impact stability, dissolution behaviour and downstream manufacturability.

“THE MECHANICAL ENERGY INTRODUCED DURING SIZE REDUCTION CAN ALTER THE SOLID-STATE FORM OF THE API BY WAY OF POLYMORPHIC TRANSITIONS, PARTIAL AMORPHISATION AND SURFACE DISORDER.”

Solid-state characterisation is therefore not a downstream confirmation step. It is the foundation for understanding how a drug substance behaves under processing conditions. When integrated early and applied consistently, it can reduce development risk, support scale-up and aid micronisation in delivering the intended performance without unintended consequences.

WHY MICRONISATION INTRODUCES RISK WITHOUT SOLID-STATE CONTROL

Micronisation is often approached as a mechanical step. In practice, it is a material-sensitive process where particle collisions, shear forces and localised temperature increases can disrupt and change crystal structure. These process-induced effects may include:

- Partial amorphisation not fully captured by bulk measurements
- Polymorphic conversion triggered by stress or humidity exposure
- Surface disorder that alters dissolution kinetics
- Recrystallisation during storage or downstream processing.

Figure 1 shows a goniometer being used for powder X-ray diffraction (PXRD), a high-precision mechanical technique that is used to accurately position and rotate the sample, X-ray source and detector to measure diffraction angles and provide an output used to confirm crystalline form.

Mechanical activation during milling has been shown to have the potential to generate disordered regions that can later recrystallise, leading to variability in product performance.^{3,4} Without a defined solid-state strategy, these risks often emerge later, during scale-up or stability studies, when mitigation is more complex and costly.




Figure 1: Goniometer used in PXRD data collection.

“MECHANICAL ACTIVATION DURING MILLING HAS BEEN SHOWN TO HAVE THE POTENTIAL TO GENERATE DISORDERED REGIONS THAT CAN LATER RECRYSTALLISE, LEADING TO VARIABILITY IN PRODUCT PERFORMANCE.”

DESIGNING A MICRONISATION STRATEGY

A successful micronisation strategy begins with the material, not the equipment. The key considerations include:

- Crystallinity, polymorphic form and thermal behaviour
- Sensitivity to moisture and mechanical stress
- Likelihood of phase transitions under process conditions
- Impact of downstream unit operations such as blending and tableting.

Materials that appear stable under ambient conditions may respond differently under high-energy milling environments. Solid-state transformations under stress are well established and can directly affect dissolution, stability and manufacturability.⁴

A risk-based approach requires linking material attributes to process parameters early in development. This includes defining critical quality attributes (CQAs), understanding process sensitivities and building data that translate at scale.

This approach aligns with quality-by-design principles, where material understanding drives process design rather than trial-and-error optimisation.⁵

MAINTAINING CONTROL OF SOLID-STATE PROPERTIES ACROSS THE PRODUCT LIFECYCLE

Solid-state risk does not disappear after micronisation. It evolves across the entire development lifecycle:

Early Development

Early development defines the foundation for success. The key questions include:

- Which polymorph, hydrate or solvate is most stable under process-relevant conditions?
- How does the material respond to thermal, mechanical and humidity stress?
- What particle size is required to achieve the target dissolution profile?

The selected solid-state form directly influences dissolution rate, stability and manufacturability.⁶ Decisions made at this stage determine the robustness of the micronisation strategy.

Analytical Technique	What It Measures	Why It Matters
X-ray Diffraction	Crystal structure, polymorph identity, amorphous content	Verifies solid-state form and crystalline structure of product
Differential Scanning Calorimetry	Thermal transitions (melting point, glass transition, recrystallisation)	Detects polymorphic changes and amorphisation introduced during processing
Thermogravimetric Analysis	Weight changes from moisture loss, decomposition or solvent evaporation	Assesses thermal stability and drying behaviour
Dynamic Vapour Sorption	Moisture or solvent uptake profiles under controlled humidity	Evaluates hygroscopicity and supports formulation and storage strategy
Specific Surface Area	Surface area per unit mass	Links surface properties to dissolution, flow and particle interactions
Particle Size Distribution	Distribution of particle sizes (e.g. via laser diffraction, dynamic light scattering)	Defines a critical quality attribute impacting dissolution and content uniformity
Scanning Electron Microscopy	Particle morphology, surface texture, agglomeration	Correlates physical structure with processing behaviour and performance

Table 1: Overview of key solid-state analytical techniques used to characterise APIs before and after micronisation.

“AS DEVELOPMENT PROGRESSES, CONSISTENCY BECOMES CRITICAL. VARIABILITY IN PARTICLE SIZE OR SOLID-STATE FORM CAN LEAD TO MEASURABLE DIFFERENCES IN DISSOLUTION AND BIOAVAILABILITY.”

Clinical Development

As development progresses, consistency becomes critical. Variability in particle size or solid-state form can lead to measurable differences in dissolution and bioavailability. Solid-state characterisation during clinical development can assist:

- Stability of the selected polymorph during micronisation
- Control of process-induced amorphous content
- Consistent particle size distribution (PSD) across batches and scales
- Reproducible material performance.

These controls are essential to maintain alignment between clinical and commercial material.⁶

Process Validation

During validation, the focus shifts from understanding to control. Solid-state characterisation can support this by:

- Confirming equivalence between input and micronised material
- Establishing acceptable process ranges
- Linking material attributes to product performance
- Supporting regulatory submissions.

Polymorphic form and solid-state properties are recognised as CQAs due to their direct impact on drug product performance.⁶

Commercial Supply

Even after process validation, solid-state risks remain. Changes in raw materials, environment or equipment can introduce variability that affects product quality. Ongoing monitoring can ensure:

- Long-term polymorphic stability
- Control of amorphous content
- Consistency in PSD and morphology
- Early detection of process drift.

Maintaining this control is essential for attaining consistent product performance and patient safety.

BUILDING A MULTI-TECHNIQUE ANALYTICAL STRATEGY

No single analytical method can fully characterise solid-state behaviour. A combination of orthogonal techniques is required to build confidence in form and stability. Table 1 provides an overview of key solid-state analytical techniques used to characterise APIs before and after micronisation.

An example of the dynamic vapour sorption (DVS) sample vessel used to determine mass changes in real time under varying humidity is shown in Figure 2.

Each technique contributes a different perspective. Together, they provide a complete understanding of how the material responds to micronisation and subsequent processing.



Figure 2: A DVS sample vessel used to determine mass changes in real time under varying humidity levels.

Applying Solid-State Thinking Beyond the API

Solid-state considerations extend beyond the API to intermediates and finished dosage forms.

API Evaluation

- Confirm polymorphic form and crystallinity
- Assess stability and moisture sensitivity
- Define PSD and morphology.

Intermediate Evaluation

- Monitor transformations during processing
- Evaluate API–excipient compatibility
- Detect amorphous content.

Final Product Evaluation

- Confirm API form within the formulation
- Verify stability over shelf life
- Link solid-state properties to dissolution and bioavailability.

This integrated approach ensures continuity in characterisation from raw materials to final product and reduces the risk of late-stage surprises.⁶

“SOLID-STATE CONSIDERATIONS EXTEND BEYOND THE API TO INTERMEDIATES AND FINISHED DOSAGE FORMS.”



America

DDF Summit

Drug Delivery & Formulation

14-15 SEPTEMBER 2026 • BOSTON

Where Pharma Innovation Meets Practical Solutions

This is your chance to join scientific leaders from across the US. With two action-packed days, and expert presentations exploring the latest in small molecules, biologics, combination products, and innovative delivery systems, it's an event you don't want to miss!

USE CODE: DDFBOSTON TO GET 20% OFF

Why attend?

- Hear from senior industry leaders about answers to the challenges you face day-to-day
- Develop contacts with like-minded individuals and enhance established connections
- Meet with leading industry solutions providers who can help you with your business

Register Now at
www.ddfsummit.com

Copyright © 2026 Furness Publishing Ltd

ONdrugDelivery • Issue 187

17

“THE UNDERLYING SOLID-STATE PROPERTIES OF THE MATERIAL ULTIMATELY DETERMINE STABILITY, PERFORMANCE AND MANUFACTURABILITY.”

CONCLUSION

Micronisation and particle engineering are powerful tools, but their success depends on more than achieving a target particle size. The underlying solid-state properties of the material ultimately determine stability, performance and manufacturability.

A development strategy that integrates solid-state characterisation before and after micronisation, supported by process understanding and lifecycle control, reduces risk and enables predictable scale-up.

For sponsors, the ability to combine particle engineering expertise with deep solid-state knowledge is critical. It enables informed decision-making, reduces development uncertainty and supports consistent product quality from early development through to commercialisation.

REFERENCES

1. Amidon GL et al, “A theoretical basis for a biopharmaceutical drug classification: The correlation of *in vitro* drug product dissolution and *in vivo* bioavailability”. *Pharm Res*, 1995, Vol 12(3), pp 413–420.
2. Csicsák D et al, “The Effect of the Particle Size Reduction on the Biorelevant Solubility and Dissolution of Poorly Soluble Drugs with Different Acid-Base Character”. *Pharmaceutics*, 2023, Vol 15(1), art 278.
3. Descamps M, Willart JF, “Perspectives on the amorphisation/milling relationship in pharmaceutical materials”. *Adv Drug Deliv Rev*, 2016, Vol 100, pp 51–66.
4. Willart JF, Descamps M, “Solid state amorphization of pharmaceuticals”. *Mol Pharm*, 2008, Vol 5(6), pp 905–920.
5. Yu LX, “Pharmaceutical quality by design: Product and process development, understanding, and control”. *Pharm Res*, 2008, Vol 25(4), pp 781–791.
6. Giron D, Mutz M, Garnier S, “Solid-state of pharmaceutical compounds”. *J Therm Anal Calorim*, 2004, Vol 77, pp 709–747.



**Dr Derek
Beauchamp**

Derek Beauchamp, PhD, Account Executive of Micronisation and Particle Engineering at Catalent, has over 20 years of pharmaceutical industry experience, with substantial expertise and knowledge in crystallography, crystallisation development and solid-state chemistry of small molecule APIs. In his role, he is responsible for scoping new opportunities, providing scientific and technical leadership for the micronisation of small molecule APIs across developmental and commercial solutions. Dr Beauchamp holds a PhD in Supramolecular Chemistry from the University of Windsor in Ontario, Canada.

T: +1 734 312 0350
E: derek.beauchamp@catalent.com



John Kusnierz

John Kusnierz is the Laboratory Manager at Catalent in Malvern (PA, US), where he leads analytical operations supporting pharmaceutical development and processing. His professional background in analytical method development and validation dates back to 2016, with experience supporting client needs across excipient, raw material, drug substance and drug product. Throughout his career, Mr Kusnierz has supported materials at all stages of the product lifecycle – from early R&D and clinical trials through process validation and commercial manufacturing. In his current role, he is responsible for analytical site leadership, programme development and the execution of all analytical work, including solid-state testing. He holds a BS degree in Marine Science from Kutztown University of Pennsylvania, US.

E: john.kusnierz@catalent.com

Catalent Pharma Services

333 Phoenixville Pike, Malvern, PA 19355, United States
www.catalent.com

WHERE CONTENT MEETS INTELLIGENCE



Championing
the missions that matter™

Tailored particle-size solutions to move your program forward

- Micronization & milling tailored to your molecule
- Scalable development-to-commercialization execution
- Proven processes for reliable yield
- Dedicated US & UK sites with expert support

MICRONIZATION & MILLING

HIGH-POTENCY & COMPLEX APIs

ANALYTICAL TESTING

©2026 Catalent, Inc. All rights reserved.

Contact us

US +1 888 SOLUTION (765-8846) EUR 00800 88 55 6178
catalent.com/micronization

Catalent
Pharma Services™